

OUR STANDARDS

- High quality services
- Excellent customer service
- Urinalysis Drug Testing performed on Immunoassay Analyzer locally.
- Screening accuracy rate of 97%-100% for rapid and reliable results.
- Comprehensive and cost effective drug and alcohol testing options.
- All urinalysis screening tests are completely processed within minutes and results are sent within 24 hours.
- Provides value-focused substance abuse screening products and lab services to governmental, criminal justice and treatment facilities.
- Committed to respecting your privacy to ensure the confidentiality of our clients.

LOCATION & CONTACT INFORMATION

Troup County Drug Screening Laboratory is located in the Troup County Government Center, on the 2nd Floor, Suite 2921, 100 Ridley Avenue, Lagrange Georgia.

Upon entrance, take the elevator to the second floor and follow the signs to the Troup County Drug Lab.



TROUP COUNTY
BOARD OF COMMISSIONERS

For more information on pricing, please contact us via phone or email.

Phone: 706-298-3615

Other: 706-298-3790

Email: mparmer@troupcountyga.gov

TROUP COUNTY DRUG SCREENING LABORATORY

Laboratory Information and Services



TROUP COUNTY
GEORGIA



Troup County Board Of Commissioners

Troup County Drug Screening Laboratory

INTRODUCTION

Established in 2011, the Troup County Drug Screening Laboratory offers urinalysis testing through Thermofisher analyzers, which uses Syva EMIT technology for high quality and reliable results. EMIT Assays are the screening tests most commonly used by labs certified by SAMHSA, and they have a screening accuracy rate of 97%-100%. Using an analyzer decreases turn around time from hours to minutes and increases drug testing efficiency.



With greater efficiencies, a sustainability and financial independence is within reach.

TESTING

The list below will provide you the EMT Assays that our lab uses, along with their respective cut offs.

- Amphetamine- 500ng/mL
- Barbiturates- 200ng/mL
- Benzodiazepine- 200ng/mL
- Cocaine- 150ng/mL
- Opiates- 300ng/mL
- Oxycodone- 100ng/mL
- THC- 50ng/mL
- ETG Metabolite (Alcohol)- 500ng/mL
- Ecstasy- 500ng/mL
- Methodon- 300ng/mL
- 6-Acelymorphine Assay- 10ng/mL
- Buprenophine- 5ng/mL
- Validity Test Creatine- 20ng/mL
- Fentanyl – 100 ng/mL
- Kratom – 50 ng/mL

ADDITIONAL SERVICES

- Hair Testing - 5 Panel
- Sweat Patch Testing- 5 Panel
- Instant Oral/ Saliva Testing- 5 Panel
- DNA Testing
- Confirmation Testing (Sent Off)- Method of testing includes: Liquid Chromatography Mass Spectrometry (LCMS) and Gas Chromatography Mass Spectrometry (GCMS)

HOW IT WORKS

- All screenings follow chain of custody procedures for sample collection and handling. All results are sent via email or fax to the requested party upon arrival to our facility.
- 5 major drug classes: Amphetamine, Cocaine, Opiates, PCP, and THC confirmed by GCMS or LCMS/MS
- In particular, the Indiko Plus analyzer uses EMIT drug-of-abuse test, also known as immunoassays. This type of test uses antibodies to detect the presence of drugs or other substances in urine. The antibodies in each assay are designed to react only with the drug being tested. If the drug being tested is present at or above the standard cut off, a chemical reaction occurs that causes a positive result.